

2020 Price list - BDE single-component solutions

BDE solutions							
CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034/ ISO 17025	PRICE, A \$
MS-BDE11ZT1	BDE 1 [CAS:7025-06-1] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	135
MS-BDE101ZT10	BDE 1 [CAS:7025-06-1] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	135
MS-BDE21ZT1	BDE 2 [CAS:6876-00-2] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	135
MS-BDE201ZT10	BDE 2 [CAS:6876-00-2] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	135
MS-BDE31ZT1	BDE 3 [CAS:101-55-3] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE301ZT10	BDE 3 [CAS:101-55-3] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE51ZT1	BDE 5 [CAS:446254-14-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	370
MS-BDE501ZT10	BDE 5 [CAS:446254-14-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	370
MS-BDE71ZT1	BDE 7 [CAS:171977-44-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	135
MS-BDE701ZT10	BDE 7 [CAS:171977-44-9] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	135
MS-BDE81ZT1	BDE 8 [CAS:147217-71-8] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE801ZT10	BDE 8 [CAS:147217-71-8] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE1001ZT10	BDE 10 [CAS:51930-04-2] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE101ZT1	BDE 10 [CAS:51930-04-2] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE121ZT1	BDE 12 [CAS:189084-59-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	75
MS-BDE1201ZT10	BDE 12 [CAS:189084-59-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	75
MS-BDE131ZT1	BDE 13 [CAS:83694-71-7] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	180
MS-BDE1301ZT10	BDE 13 [CAS:83694-71-7] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	190
MS-BDE141ZT1	BDE 14 [CAS:46438-88-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	375
MS-BDE1401ZT10	BDE 14 [CAS:46438-88-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	375
MS-BDE151ZT1	BDE 15 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE1501ZT10	BDE 15 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE171ZT1	BDE 17 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE1701ZT10	BDE 17 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE181ZT1	BDE 18 [CAS:407606-55-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	335
MS-BDE1801ZT10	BDE 18 [CAS:407606-55-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	335
MS-BDE191ZT1	BDE 19 [CAS:147217-73-0] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	175
MS-BDE1901ZT10	BDE 19 [CAS:147217-73-0] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	175
MS-BDE211ZT1	BDE 21 [CAS:N/A] 100mg/l in Iso-octane	1	ml	24	BDE solution	Yes	375
MS-BDE2101ZT10	BDE 21 [CAS:N/A] 10mg/l in Iso-octane	10	ml	24	BDE solution	Yes	375
MS-BDE221ZT1	BDE 22 [CAS:446254-15-5] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	375
MS-BDE2201ZT10	BDE 22 [CAS:446254-15-5] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	375
MS-BDE231ZT1	BDE 23 [CAS:446254-16-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	375
MS-BDE2301ZT10	BDE 23 [CAS:446254-16-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	375
MS-BDE241ZT1	BDE 24 [CAS:NA] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	375
MS-BDE2401ZT10	BDE 24 [CAS:NA] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	375

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034/ ISO 17025	PRICE, A \$
MS-BDE251ZT1	BDE 25 [CAS:147217-77-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE2501ZT10	BDE 25 [CAS:147217-77-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE261ZT1	BDE 26 [CAS:337513-75-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	290
MS-BDE2601ZT10	BDE 26 [CAS:337513-75-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	290
MS-BDE271ZT1	BDE 27 [CAS:337513-53-8] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	175
MS-BDE2701ZT10	BDE 27 [CAS:337513-53-8] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	175
MS-BDE281ZT1	BDE 28 [CAS:41318-75-6] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE2801ZT10	BDE 28 [CAS:41318-75-6] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	135
MS-BDE291ZT1	BDE 29 [CAS:337513-56-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	370
MS-BDE2901ZT10	BDE 29 [CAS:337513-56-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	370
MS-BDE301ZT1	BDE 30 [CAS:155999-95-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	175
MS-BDE3001ZT10	BDE 30 [CAS:155999-95-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	175
MS-BDE311ZT1	BDE 31 [CAS:65075-08-3] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	340
MS-BDE3101ZT10	BDE 31 [CAS:65075-08-3] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	335
MS-BDE321ZT1	BDE 32 [CAS:189084-60-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	145
MS-BDE3201ZT10	BDE 32 [CAS:189084-60-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	145
MS-BDE331ZT1	BDE 33 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	246
MS-BDE3301ZT10	BDE 33 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	246
MS-BDE341ZT1	BDE 34 [CAS:446254-17-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE3401ZT10	BDE 34 [CAS:446254-17-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE351ZT1	BDE 35 [CAS:147217-80-9] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	134
MS-BDE3501ZT10	BDE 35 [CAS:147217-80-9] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	134
MS-BDE371ZT1	BDE 37 [CAS:147217-81-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE3701ZT10	BDE 37 [CAS:147217-81-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE381ZT1	BDE 38 [CAS:337513-54-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE3801ZT10	BDE 38 [CAS:337513-54-9] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE391ZT1	BDE 39 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE3901ZT10	BDE 39 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE411ZT1	BDE 41 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE4101ZT10	BDE 41 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE421ZT1	BDE 42 [CAS:446254-18-8] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE4201ZT10	BDE 42 [CAS:446254-18-8] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE441ZT1	BDE 44 [CAS:446254-20-2] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE4401ZT10	BDE 44 [CAS:446254-20-2] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE461ZT1	BDE 46 [CAS:446254-22-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	361
MS-BDE4601ZT10	BDE 46 [CAS:446254-22-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	361
MS-BDE471ZT1	BDE 47 [CAS:5436-43-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	135
MS-BDE4701ZT10	BDE 47 [CAS:5436-43-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	70
MS-BDE481ZT1	BDE 48 [CAS:337513-55-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	350

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034\ ISO 17025	PRICE, A \$
MS-BDE4801ZT10	BDE 48 [CAS:337513-55-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	350
MS-BDE491ZT1	BDE 49 [CAS:243982-82-3] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	310
MS-BDE4901ZT10	BDE 49 [CAS:243982-82-3] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	310
MS-BDE501ZT1	BDE 50 [CAS:446254-23-5] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	180
MS-BDE5001ZT10	BDE 50 [CAS:446254-23-5] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	180
MS-BDE511ZT1	BDE 51 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	180
MS-BDE5101ZT10	BDE 51 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	180
MS-BDE521ZT1	BDE 52 [CAS:446254-24-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE5201ZT10	BDE 52 [CAS:446254-24-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE531ZT1	BDE 53 [CAS:446254-25-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE5301ZT10	BDE 53 [CAS:446254-25-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE541ZT1	BDE 54 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE5401ZT10	BDE 54 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE551ZT1	BDE 55 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE5501ZT10	BDE 55 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE561ZT1	BDE 56 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE5601ZT10	BDE 56 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE571ZT1	BDE 57 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE5701ZT10	BDE 57 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE581ZT1	BDE 58 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE5801ZT10	BDE 58 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE601ZT1	BDE 60 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	362
MS-BDE6001ZT10	BDE 60 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	362
MS-BDE611ZT1	BDE 61 [CAS:446254-32-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	362
MS-BDE6101ZT10	BDE 61 [CAS:446254-32-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	362
MS-BDE62ZT1	BDE 62 [CAS:446254-33-7] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	290
MS-BDE6201ZT10	BDE 62 [CAS:446254-33-7] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	290
MS-BDE631ZT1	BDE 63 [CAS:446254-34-8] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE6301ZT10	BDE 63 [CAS:446254-34-8] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE641ZT1	BDE 64 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE6401ZT10	BDE 64 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE651ZT1	BDE 65 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE6501ZT10	BDE 65 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE661ZT1	BDE 66 [CAS:189084-61-5] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	185
MS-BDE6601ZT10	BDE 66 [CAS:189084-61-5] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	185
MS-BDE671ZT1	BDE 67 [CAS:446254-37-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	352
MS-BDE6701ZT10	BDE 67 [CAS:446254-37-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	352
MS-BDE681ZT1	BDE 68 [CAS:446254-38-2] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	310
MS-BDE6801ZT10	BDE 68 [CAS:446254-38-2] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	310

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034/ ISO 17025	PRICE, A \$
MS-BDE691ZT1	BDE 69 [CAS:327185-09-1] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	185
MS-BDE6901ZT10	BDE 69 [CAS:327185-09-1] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	185
MS-BDE701ZT1	BDE 70 [CAS:446254-39-3] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE7001ZT10	BDE 70 [CAS:446254-39-3] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	440
MS-BDE711ZT1	BDE 71 [CAS:189084-62-6] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	184
MS-BDE7101ZT10	BDE 71 [CAS:189084-62-6] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	180
MS-BDE731ZT1	BDE 73 [CAS:446254-41-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE7301ZT10	BDE 73 [CAS:446254-41-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE741ZT1	BDE 74 [CAS:446254-42-8] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE7401ZT10	BDE 74 [CAS:446254-42-8] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE751ZT1	BDE 75 [CAS:189084-63-7] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	184
MS-BDE7501ZT10	BDE 75 [CAS:189084-63-7] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	180
MS-BDE761ZT1	BDE 76 [CAS:446254-43-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE771ZT1	BDE 77 [CAS:93703-48-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	184
MS-BDE7701ZT10	BDE 77 [CAS:93703-48-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	184
MS-BDE781ZT1	BDE 78 [CAS:446254-45-1] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE7801ZT10	BDE 78 [CAS:446254-45-1] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE791ZT1	BDE 79 [CAS:446254-48-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	310
MS-BDE7901ZT10	BDE 79 [CAS:446254-48-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	310
MS-BDE801ZT1	BDE 80 [CAS:103173-66-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE8001ZT10	BDE 80 [CAS:103173-66-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE811ZT1	BDE 81 [CAS:446254-50-8] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE8101ZT10	BDE 81 [CAS:446254-50-8] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE821ZT1	BDE 82 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE8201ZT10	BDE 82 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE831ZT1	BDE 83 [CAS:446254-51-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE8301ZT10	BDE 83 [CAS:446254-51-9] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE851ZT1	BDE 85 [CAS:182346-21-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE8501ZT10	BDE 85 [CAS:182346-21-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE861ZT1	BDE 86 [CAS:446254-53-1] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE8601ZT10	BDE 86 [CAS:446254-53-1] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE871ZT1	BDE 87 [CAS:446254-54-2] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE8701ZT10	BDE 87 [CAS:446254-54-2] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	440
MS-BDE881ZT1	BDE 88 [CAS:446254-55-3] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	320
MS-BDE8801ZT10	BDE 88 [CAS:446254-55-3] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	320
MS-BDE891ZT1	BDE 89 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE8901ZT10	BDE 89 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	440
MS-BDE901ZT1	BDE 90 [CAS:446254-57-5] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9001ZT10	BDE 90 [CAS:446254-57-5] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034 ISO 17025	PRICE, A \$
MS-BDE911ZT1	BDE 91 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9101ZT10	BDE 91 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE921ZT1	BDE 92 [CAS:446254-59-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9201ZT10	BDE 92 [CAS:446254-59-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE931ZT1	BDE 93 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	445
MS-BDE9301ZT10	BDE 93 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	445
MS-BDE941ZT1	BDE 94 [CAS:446254-61-1] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9401ZT10	BDE 94 [CAS:446254-61-1] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE951ZT1	BDE 95 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9501ZT10	BDE 95 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE961ZT1	BDE 96 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9601ZT10	BDE 96 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE971ZT1	BDE 97 [CAS:446254-64-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9701ZT10	BDE 97 [CAS:446254-64-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE981ZT1	BDE 98 [CAS:38463-82-0] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE9801ZT10	BDE 98 [CAS:38463-82-0] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE991ZT1	BDE 99 [CAS:60348-60-9] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE9901ZT10	BDE 99 [CAS:60348-60-9] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1001ZT1	BDE 100 [CAS:189084-64-8] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE10001ZT10	BDE 100 [CAS:189084-64-8] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1011ZT1	BDE 101 [CAS:446254-65-5] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE10101ZT10	BDE 101 [CAS:446254-65-5] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1021ZT1	BDE 102 [CAS:446254-66-6] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE10201ZT10	BDE 102 [CAS:446254-66-6] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1031ZT1	BDE 103 [CAS:446254-67-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	285
MS-BDE10301ZT10	BDE 103 [CAS:446254-67-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	285
MS-BDE1041ZT1	BDE 104 [CAS:446254-68-8] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE10401ZT10	BDE 104 [CAS:446254-68-8] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1051ZT1	BDE 105 [CAS:373594-78-6] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE10501ZT10	BDE 105 [CAS:373594-78-6] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1061ZT1	BDE 106 [CAS:446254-69-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE10601ZT10	BDE 106 [CAS:446254-69-9] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1071ZT1	BDE 107 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE10701ZT10	BDE 107 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1081ZT1	BDE 108 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	320
MS-BDE10801ZT10	BDE 108 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	320
MS-BDE1091ZT1	BDE 109 [CAS:446254-72-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	445
MS-BDE10901ZT10	BDE 109 [CAS:446254-72-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	445
MS-BDE1101ZT1	BDE 110 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034 ISO 17025	PRICE, A \$
MS-BDE11001ZT10	BDE 110 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1111ZT1	BDE 111 [CAS:446254-74-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE11101ZT10	BDE 111 [CAS:446254-74-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1121ZT1	BDE 112 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE11201ZT10	BDE 112 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1131ZT1	BDE 113 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE11301ZT10	BDE 113 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1141ZT1	BDE 114 [CAS:446254-77-9] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE11401ZT10	BDE 114 [CAS:446254-77-9] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1151ZT1	BDE 115 [CAS:446254-78-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	320
MS-BDE11501ZT10	BDE 115 [CAS:446254-78-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	320
MS-BDE1161ZT1	BDE 116 [CAS:189084-65-9] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE11601ZT10	BDE 116 [CAS:189084-65-9] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1171ZT1	BDE 117 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE11701ZT10	BDE 117 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1181ZT1	BDE 118 [CAS:446254-80-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE11801ZT10	BDE 118 [CAS:446254-80-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1191ZT1	BDE 119 [CAS:189084-66-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE11901ZT10	BDE 119 [CAS:189084-66-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1201ZT1	BDE 120 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	320
MS-BDE12001ZT10	BDE 120 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	320
MS-BDE1211ZT1	BDE 121 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE12101ZT10	BDE 121 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE1221ZT1	BDE 122 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE12201ZT10	BDE 122 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1231ZT1	BDE 123 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE12301ZT10	BDE 123 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1241ZT1	BDE 124 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE12401ZT10	BDE 124 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1251ZT1	BDE 125 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE12501ZT10	BDE 125 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1261ZT1	BDE 126 [CAS:366791-32-4] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE12601ZT10	BDE 126 [CAS:366791-32-4] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1271ZT1	BDE 127 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE12701ZT10	BDE 127 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1281ZT1	BDE 128 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE12801ZT10	BDE 128 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1291ZT1	BDE 129 [CAS:NA] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE12901ZT10	BDE 129 [CAS:NA] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034 ISO 17025	PRICE, A \$
MS-BDE15301ZT10	BDE 153 [CAS:68631-49-2] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1541ZT1	BDE 154 [CAS:207122-15-4] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE15401ZT10	BDE 154 [CAS:207122-15-4] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1551ZT1	BDE 155 [CAS:35854-94-5] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	187
MS-BDE15501ZT10	BDE 155 [CAS:35854-94-5] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	187
MS-BDE1561ZT1	BDE 156 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE15601ZT10	BDE 156 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1571ZT1	BDE 157 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE15701ZT10	BDE 157 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1581ZT1	BDE 158 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE15801ZT10	BDE 158 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1591ZT1	BDE 159 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE15901ZT10	BDE 159 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1611ZT1	BDE 161 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE16101ZT10	BDE 161 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1621ZT1	BDE 162 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE16201ZT10	BDE 162 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1631ZT1	BDE 163 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE16301ZT10	BDE 163 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1641ZT1	BDE 164 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE16401ZT10	BDE 164 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1651ZT1	BDE 165 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE16501ZT10	BDE 165 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1661ZT1	BDE 166 [CAS:189084-58-0] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	185
MS-BDE16601ZT10	BDE 166 [CAS:189084-58-0] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	185
MS-BDE1671ZT1	BDE 167 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE16701ZT10	BDE 167 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE1681ZT1	BDE 168 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE16801ZT10	BDE 168 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1691ZT1	BDE 169 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE16901ZT10	BDE 169 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE1701ZT1	BDE 170 [CAS:327185-13-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE17001ZT10	BDE 170 [CAS:327185-13-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE1711ZT1	BDE 171 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE17101ZT10	BDE 171 [CAS:N/A] 10ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE1721ZT1	BDE 172 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE17201ZT10	BDE 172 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE1741ZT1	BDE 174 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	490
MS-BDE17401ZT10	BDE 174 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	490

CODE	DESCRIPTION	VOL	UNITS	SHELF LIFE	CLASS	ISO 17034 ISO 17025	PRICE, A \$
MS-BDE1951ZT1	BDE 195 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	440
MS-BDE1961ZT1	BDE 196 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE19601ZT10	BDE 196 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE1971ZT1	BDE 197 [CAS:117964-21-3] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	365
MS-BDE19701ZT10	BDE 197 [CAS:117964-21-3] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	365
MS-BDE1991ZT1	BDE 199 [CAS:NA] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	490
MS-BDE19901ZT10	BDE 199 [CAS:NA] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	490
MS-BDE2011ZT1	BDE 201 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE20101ZT10	BDE 201 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE2021ZT1	BDE 202 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE20201ZT10	BDE 202 [CAS:N/A] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE2031ZT1	BDE 203 [CAS:337513-72-1] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE20301ZT10	BDE 203 [CAS:337513-72-1] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE204ZT1	BDE 204 [CAS:N/A] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	360
MS-BDE20401ZT10	BDE 204 [CAS:N/A] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	360
MS-BDE2051ZT1	BDE 205 [CAS:446255-56-7] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	360
MS-BDE20501ZT10	BDE 205 [CAS:446255-56-7] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	360
MS-BDE2061ZT1	BDE 206 [CAS:63387-28-0] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE20601ZT10	BDE 206 [CAS:63387-28-0] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE2071ZT1	BDE 207 [CAS:437701-79-6] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE20701ZT10	BDE 207 [CAS:437701-79-6] 10ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE2081ZT1	BDE 208 [CAS:N/A] 100ug/ml in Iso-octane Storage conditions: Stored under normal laboratory conditions	1	ml	24	BDE solution	Yes	440
MS-BDE20801ZT10	BDE 208 [CAS:N/A] 10mg/l in Iso-octane Storage conditions: Stored under normal laboratory conditions	10	ml	24	BDE solution	Yes	440
MS-BDE2091ZT1	BDE 209 [CAS:1163-19-5] 100ug/ml in Iso-octane	1	ml	24	BDE solution	Yes	68
MS-BDE20901ZT10	BDE 209 [CAS:1163-19-5] 10ug/ml in Iso-octane	10	ml	24	BDE solution	Yes	68