

2021		Moolap Scientific P/L	Price
MOOSCI P/N°	OEM P/No.	Parts List for Spectro EOP ICP	A \$
2050-0710NT		'Tyfoon' high-TDS concentric glass nebulizer, with NebLink sample tube	690
SILQ-50-A0.5		Meinhard Ultra-Hi-Purity QUARTZ concentric nebulizer, PFA argon and sample lines, low-uptake	990
TR-30-A1	TR-30-A1	Meinhard Type A Concentric Glass Nebulizer, nom. uptake 1mL/min, 1 L/min @ 30psi	470
TR-30-C1	76060510	Meinhard Type C Concentric glass nebulizer, nom. uptake 1mL/min, 1 L/min @ 30psi	640
TR-30-K2	TR-30-K2	Meinhard Concentric glass nebulizer, Type K, nom. uptake 2mL/min, 1 L/min @ 30psi	540
EP-2050-0700N	76060510	C-Type Concentric Glass Nebulizer, with 'NebLink' sample tube	540
EP-8100-2051		'PEARL' Cleaning Device for <i>all types</i> of concentric glass and QUARTZ nebulizers	290
EP-2045-0700		'PERGO' Argon Humidifier, with tubing and connectors	1480
3150-0720	PMM 4000	Burgener PEEK <i>MiraMist</i> Parallel-Path Nebulizer for Spectro ICP	1190
EP-2050-4251		'Boreal' Medium-Uptake 'No-Break' concentric polymer nebulizer	870
2550-2024	ME2040	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	1370
2030-2028	ME2030	PFA-HS HF-res. High-Solids concentric neb, with 400uL/min detachable sample uptake tube	1490
2550-2028	ME4040	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	720
EP-2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	110
EP-2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	110
EP-2050-092S	75060502	Fixed Cross-Flow for Spectro, 34mm Inner Diameter	830
EP-2051-020S	48205036	Spectro Modified Lichte MDSN Nebulizer	830
ML175026	48105001	Cyclonic spray chamber, top socket, pump drain	295
ML175026ES		Cyclonic spray chamber, top socket, pump drain with Easy Seal nebulizer inlet and PFA cap	535
ML175026EOP		Spectro Cyclonic Spray Chamber, end-on plasma (EOP), Pump Drain	355
ML175026EOP-ES		Spectro Cyclonic Spray Chamber, Top Socket, Pump Drain with Easy-Seal	535
ML175026Q		Quartz Cyclonic Spray Chamber, Socket	575
ML175026QJ	USA00850	Quartz Cyclonic Spray Chamber, Jacketed, Threaded Ports	885
ME31671011	ES-3167-1011	Quartz Cyclonic Spray Chamber for PC 3 Peltier cooler	585
EP-2075-0120	48105003	Spectro Small Conical Spray Chamber, with Impact Sphere, vertical exitm ML175012	215
EP-2075-0130	48105051	Spectro Glass Spray Chamber, Double Pass 90° Exit, ML175019	295
EP-3167-1001	ES-3167-1001	Quartz cyclonic spray chamber, can be used with PC3 Peltier cooling unit	995
EP-2075-0140	48105006	Spectro Large Double Pass Spray Chamber, RS12/5 outlet, ML175014	290
EP-2256-5470	ME22565470	47-mm PFA spray chamber for Spectro ICP-MS, for use with ESI torch	1475
EP-4167-1011	ME41671011	PC3 Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply	3890
EP-2250-0470	ES-2250-0470	PFA Sapphire HF-resistant Kit for Spectro CIROS, with PFA-ST nebulizer	3990
EP-2075-0170	48105062	Spectro Elbow Adapter, EOP torch to Chamber	70
EP-2075-0180	48105002	Spectro Elbow Adapter Slurry, EOP Torch to Chamber	70
EP-2075-0300		Spectro Drain/Pressure Lock	245
EP-2075-1702	ES-1702-0000	Quartz Torch for demountable <i>axial-view</i> CIROS	490
EP-2075-0210	48105008	Spectro Quartz Injector/Aerosol Tube 1.5mm, ML175021	160
EP-2075-0220	48105010	Spectro Quartz Injector/Aerosol Tube 1.8mm, ML175022	160
EP-2075-0230	48105009	Spectro Quartz Injector/Aerosol Tube 2.5mm, ML175023	160
EP-1980-0180	ES-1980-0180	Sapphire Injector in PFA adaptor, ID = 1,8 mm	1090
EP-1980-0300	ES-1980-0300	Sapphire Injector in PFA adaptor, ID = 3 mm	1090
EP-1981-0200	ES-1981-0200	Platinum Injector in PFA adaptor, ID = 2 mm	2570
EP-2075-0400	48105052	Spectro EOP Flared End Quartz Torch 2.5mm injector, 4mm side arms, ML175040	520
EP-2075-0400F	75160526	Spectro EOP Flared End Quartz Torch 2.5mm injector, with fittings	580
ML175041		Spectro Quartz Torch, end-on plasma (EOP), Axial, Fixed Flared End, ID=2.2mm, ML175041	520
EP-2075-0410	48105052A	Spectro EOP Fixed Flared End Quartz Torch, 1.8mm ID injector, ML175042	520
EP-2075-0420	48206007	Spectro EOP Demountable Torch Tube Assembly, L=80mm, ML175044	340
ML175047	75260533	Quartz torch, end-on plasma (EOP), Axial, with sheath gas, fi xed 3.2 mm injector	695
ML175048	75160537	Quartz torch, end-on plasma (EOP), Axial, with sheath gas, fi xed 2.5 mm injector	695
ML175049	75160549	Quartz torch, end-on plasma (EOP), Axial, with sheath gas, fi xed 1.8 mm injector	695
ME1702EOP		Quartz torch, end-on plasma (EOP), Axial, demountable with flare, L=80mm	430
EP-2075-0440		Spectro EOP Demountable Flared End Torch	370
EP-2075-0460	48105058	Quartz Injector Tube, 2.5mm ID	160
EP-2075-0480	75160537	Spectro EOP Fixed Sheath Gas System Torch, OEM 75160549, ML175049	695

EP-2075-0550		Gas Connectors with Part Nos. 44302072 and 443302077	55
MOOSCI P/N°		Autosampler accessories	A \$
EP-5016-3250		Spectro AS400 Autosampler Probe, 0.25-mm ID, ~ 100 ul/min.	390
EP-5017-0750-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,8 mm	390
EP-5017-0500-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,5 mm	390
MOOSCI P/N°	OEM P/No.	Adapters and peri-pump tubing	A \$
ERT-040		EPOND 'RubeTubes', Tygon-core sample tube connectors, packet of 40	165
MM175055		Tubing and gas connectors (44302072 & 44302077)	38
Tube-PVCSD-15	47201037	2-stop Tygon pump tubing, 0.38 mm ID, 151 mm bridge distance (orange/green), 12 pc	36
Tube-PVCSD-35	47101003	2-stop Tygon pump tubing, 1.75 mm ID, 151 mm bridge distance (orange/orange), 12 pc	36
Tube-PVCSD-51	47501010	2-stop Tygon pump tubing, 1.30 mm ID, 151 mm bridge distance (grey/grey), 12 pc	36
Tube-PVCSD-69	47101005	2-stop Tygon pump tubing, 1.75 mm ID, 151 mm bridge distance (blue/green), 12 pc	36
Tube-PVCSD-15-F	4720137	2-stop Tygon pump tubing, 0.38 mm ID, 2.0 mm flared (orange/green), 12 pc	59
Tube-PVCSD-35-F	47101003	2-stop Tygon pump tubing, 0.89 mm ID, 2.0 mm flared, 151 mm bridge, (orange/orange), 12 pc	59
Tube-PVCSD-51-F	47501010	2-stop Tygon pump tubing, 1.30 mm ID, 2.0 mm flared, 151 mm bridge, (grey/grey), 12 pc	59
		<i>... custom standard solutions supplied - according to your specifications ...</i>	
www.moosci.com			